

ACTINO GEL

DENTAL ETCHING GEL

PREVESTDenPro®
THE FUTURE OF DENTISTRY

Product Description : Actino gel is a blue tinted viscous gel of 37% Phosphoric acid for etching tooth enamel and dentin conditioning or composite surfaces prior to the bonding of any composites, sealants or bonding adhesives.

Composition : Phosphoric Acid, Glycerine, Solvent, Purified Water, Fillers and Coloring Pigments.

Intended Purpose / Intended Use : It is intended to be used as enamel and dentin etching gel.

Indications for Use : Used for conditioning both dental enamel and dentine surfaces prior to the placement of an adhesive restoration.

Instructions for Use : Clean the area to be etched, rinse and dry thoroughly. Place the tip against the surface to be etched and slowly dispense the gel on the surface to lay a fine line or paint a thin layer to the areas to be etched. Allow the gel to stay for 10-15 seconds and rinse thoroughly with water spray. Dry etched tooth surface with clean, dry air. The dispensed material shall not be placed back in the container and the dose once applied or used shall not be reused.

Warnings : Avoid prolonged or repeated contact with skin (allergic contact dermatitis), and eyes. Avoid prolonged inhalation. Do not take internally. Contact with the skin may cause irritation, burns of hypersensitivity. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention. If the product comes into contact with the eyes, rinse thoroughly with fresh water and consult a physician. It should not be used with patients who have a history of severe allergic reaction to any of the components.

Precautions : The product contains phosphoric acid. It can burn and cause irritation to the tissues and therefore confine the use to areas to be etched. Fix the applicator tip to the syringe and check the flow of gel before applying on the etching surface. Avoid contact with eyes and prolonged contact with skin.

Contraindications : Its use is contraindicated in patients who are allergic to or sensitive to phosphoric acid.

Storage Conditions :

Store at temperature between 10°C(50°F) to 24°C(75°F). Close the cap carefully after use. Keep away from moisture. Shelf life not affected after opening if stored in a cool dry place.

Shelf Life : 3 years from the date of manufacturing.

Presentation :

Trial Pack : 1 x 2ml Syringe, 5 Applicator Tips (Ref : 10006-1), **Trial Pack :** 1 x 3ml Syringe, 10 Applicator Tips (Ref : 10006-2), **Trial Pack :** 1 x 3ml Syringe, 10 Applicator Tips in Box Pack (Ref : 10006-3), **Intro Pack :** 1 x 5ml Syringe, 10 Applicator Tips. (Ref : 10006), **Economy Pack :** 3 x 3ml Syringe, 30 Applicator Tips (Ref : 10006-5), **Economy Pack :** 4 x 5ml Syringe, 40 Applicator Tips (Ref : 10007), **Jumbo Pack :** 1 x 50ml Syringe Actino Gel, 2 x Empty Syringes 5 ml, 1 x 20 Applicator Tips (Ref : 10008)

Disposal conditions : To be disposed off based on the prevailing rules of the country in which it is used.

Prevest DenPro Limited

Unit II, Export Promotion Industrial Park,
(EPIP), Bari Brahmna, Jammu-181133,
India. www.prestdenpro.com
Customer Care No. : +91 8899074151
Email : customercare@prevestdenpro.com
Made in India

EC REP
Obelis s.a
Boulevard General Wahnis 53
1030 Brussels, Belgium
Tel: +32(2) 732 59 54. Fax: +32(2) 732 60 03
E-mail: mail@obelis.net

UK REP
Obelis UK Ltd
Sandford Gate, East Point
Business Park, Oxford, OX4 6LB,
UK. Phone: +44.1491.378012
E-mail: info@obelis.co.uk

CE
0123

Issue/Revision : 01/01
Date : 19-09-2025

UDI Symbol for "UNIQUE DEVICE IDENTIFIER"

LOT Symbol for "BATCH CODE" USE-BY DATE Symbol for "USE-BY DATE" NON STERILE Symbol for "NON STERILE" KEEP DRY Symbol for "KEEP DRY"

CAUTION Symbol for "CAUTION" CONSULT INSTRUCTION FOR USE Symbol for "CONSULT INSTRUCTION FOR USE" REF Symbol for "CATALOGUE NUMBER"

EC/REP Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"

DO NOT USE IF PACKAGE IS DAMAGE Symbol for "DO NOT USE IF PACKAGE IS DAMAGE" KEEP OUT OF THE REACH OF CHILDREN Symbol for "KEEP OUT OF THE REACH OF CHILDREN" MANUFACTURER Symbol for "MANUFACTURER"

TEMPERATURE LIMIT Symbol for "TEMPERATURE LIMIT" KEEP AWAY FROM SUNLIGHT Symbol for "KEEP AWAY FROM SUNLIGHT" MD Symbol for "MEDICAL DEVICE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.