

MaxGenol

Liquid for temporary cementation with zinc oxide

Product Description :

Maxgenol has an obtundent effect and reacts with zinc oxide to form a temporary cement for sedative dental dressing or lining under non-composite resin restoration.

Composition : Eugenol and Vegetable Oils.

Intended Purpose / Intended Use :

It is used as an obtundent and a temporary restorative material in combination with zinc oxide.

Indications for Use :

For temporary relief of localized dental pain. Recommended to be combined with Zinc Oxide to form a temporary sedative dressing.

Instructions for Use :

Place a drop on a cotton bud and gently apply to the inflamed area. Mix with Zinc Oxide to form a thick paste and place into a tooth cavity either as a dressing on a painful tooth or as a lining under a non composite resin restoration. The dispensed material shall not be placed back in container and the dose once applied and used shall not be reused.

Warnings :

Contact with the skin may cause irritation, burns of hypersensitivity. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue the use and seek medical attention. If the product comes into contact with the eyes, rinse thoroughly with fresh water and consult a physician. It should not be used with patients who have a history of severe allergic reaction to any of the components.

Precautions :

Avoid contact with eyes. Inhibits the polymerization of composite resin and should not be used as a dressing prior to placing a composite resin restoration.

Contraindications :

It is contraindicated for use in patients who have a history of severe allergic reaction to any of its components.

Storage Conditions :

Store in a cool and dry environment between 10°C (50°F) to 24°C (75°F). Always close the cap immediately after use. Keep away from moisture. Shelf Life is not affected after opening if stored in a cool dry place. If not to be used for a longer duration, keep in refrigerator. Allow material to reach room temperature prior to use.

Shelf Life : 3 years from the date of manufacturing.

Presentation :

Ref: 40038

1x 15ml Glass Bottle, 1 x Dispensing Dropper

Ref: 40039

1x 30ml Glass Bottle, 1 x Dispensing Dropper

Ref. : 40040

1x 110ml Glass Bottle, 1 x Dispensing Dropper

Disposal Conditions : To be disposed off based on the prevailing rules of the country in which it is used.

Prevest DenPro Limited

Unit II, Export Promotion Industrial Park,
(EPIP), Bari Brahmana, Jammu-181133,
India. www.prevestdenpro.com

Customer Care No. : + 91 8899074151

Email : customercare@prevestdenpro.com

Made in India

| | |
|-------------------|------------|
| Revision : | 04 |
| Date : | 28-09-2023 |



The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaint and reporting of unexpected events contact customer care.