

# MICRON DENTIN CONDITIONER

## POLYACRYLIC ACID CAVITY CONDITIONER

**PREVESTDenPro®**  
THE FUTURE OF DENTISTRY

**Product Description :** Micron Dentin Conditioner is a mild polyacrylic acid solution designed to remove the dentinal smear layer and to condition dentine, thus enhancing the bond between glass ionomer cement and the dentine.

**Composition :** 10% Polyacrylic Acid Solution.

**Intended Purpose / Intended Use :** It is intended to be used as pretreatment solution for dentin conditioner.

**Indications for Use :** It is used for preparing dentinal surface prior to bonding or cementing or filling.

### Instructions for Use :

1. After tooth preparation, apply Micron Dentin Conditioner to the bonding surfaces for 20 seconds using a cotton pellet or sponge.
2. Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. **DO NOT DESICCATE.** Best results are obtained when prepared surfaces appear moist (glistening).
3. Apply Micron Superior or Micron Luting or other glass ionomer as a bonding agent, base/liner or restorative material.
4. Close the bottle immediately after use.
5. The dispensed material shall not be placed back in container and the dose once applied and used shall not be reused.

**Warnings:** Avoid prolonged or repeated contact with skin (allergic contact dermatitis), and eyes. Avoid prolonged inhalation. Do not take internally. Contact with the skin may cause irritation, burns of hypersensitivity. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention. If the product comes into contact with the eyes, rinse thoroughly with fresh water and consult a physician. It should not be used with patients who have a history of severe allergic reaction to any of the components.

### Precautions :

1. In case of contact with oral tissue or skin, flush immediately with water.
  2. In case of contact with eyes, flush immediately with water and seek medical attention.
- Contraindications :** In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**Storage Conditions :** Store at temperature between 10°C(50°F) to 24°C(75°F). Close the cap carefully after use. Keep away from moisture. Shelf life not affected after opening if stored in a cool dry place.

**Shelf life :** 3 years from the date of manufacture.

### Presentation :

Ref. 30010-1

1x10 ml Dropper Bottle.

**Disposal conditions :** To be disposed off based on the prevailing rules of the country in which it is used.

**Prevest DenPro Limited**  
Unit II, EPIP, Bari Brahmna,  
Jammu-181133, India.  
www.prevestdenpro.com  
Customer Care No. : + 91 8899074151  
Email : customercare@prevestdenpro.com  
Made in India

EC REP  
Obelis s.a  
Boulevard General Wahis 53  
1030 Brussels, Belgium  
Tel: +32) 2 732 59 54, Fax: +32) 2 732 60 03  
E-mail : mail@obelis.net

CE  
0123

Revision :	06
Date :	24-01-2023

UDI Symbol for "UNIQUE DEVICE IDENTIFIER"

**LOT** Symbol for "BATCH CODE" **USE-BY DATE** Symbol for "USE-BY DATE" **NON STERILE** Symbol for "NON STERILE" **KEEP DRY** Symbol for "KEEP DRY"

**CAUTION** Symbol for "CAUTION" **CONSULT INSTRUCTION FOR USE** Symbol for "CONSULT INSTRUCTION FOR USE" **REF** Symbol for "CATALOGUE NUMBER"

**EC REP** Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"

**DO NOT USE IF PACKAGE IS DAMAGED** Symbol for "DO NOT USE IF PACKAGE IS DAMAGED" **KEEP OUT OF THE REACH OF CHILDREN** Symbol for "KEEP OUT OF THE REACH OF CHILDREN" **MANUFACTURER** Symbol for "MANUFACTURER"

**TEMPERATURE LIMIT** Symbol for "TEMPERATURE LIMIT" **KEEP AWAY FROM SUNLIGHT** Symbol for "KEEP AWAY FROM SUNLIGHT" **MD** Symbol for "MEDICAL DEVICE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.