

IFU Size : 161 H x 104 W

CeraEtch

Porcelain Etching Gel (5% Buffered Hydrofluoric Acid)



Product Description : CeraEtch is 5% Buffered Hydrofluoric acid in semi gel form. It improves the mechanical bond between porcelain and bonding agent by increasing the surface roughness.

Composition : Hydrofluoric Acid, Ammonium Fluoride, Purified Water and Thickening Agent.

Intended purpose / Intended Use : It is used as an auxiliary material in dental bonding procedure.

Indication for Use :

1. Repair of porcelain restorations.
2. Cementation of porcelain restorations.

Instruction for Use (for cementation of ceramic restoration) :

1. Protective eye wear and gloves must be worn by the practitioner.
2. Remove cap from the CeraEtch syringe and attach dispensing tip and check the etchant flow through the tip before applying. Apply etchant to the dry porcelain for 90 seconds.
3. Rinse with a copious amount of water and air dry. The etched surface should dull and frosty.
4. Apply 1– 2 thin coats of Silane-X with a brush tip applicator to the etched surface and allows to dry for 30 seconds.
5. Apply Renew Universal / Total Etch Bonding adhesive and follow the instructions for bonding and restorative with composites or cements.

Instruction for Use (for repair) :

1. Protective eye wear must be worn by the practitioner and patient and protecting gloves must be worn by the practitioner.
2. Isolate the area to be etched with a rubber dam. Apply gingival barrier around the area to be etched.
3. Remove cap from the CeraEtch syringe and attach dispensing tip and check the etchant flow through the tip before applying intra orally. Apply etchant to the dry porcelain for 90 seconds. Suction the porcelain etchant and rinse with a copious amount of water and air dry. The etched surface should dull and frosty.
4. Apply 1– 2 thin coats of Silane-X with a brush tip applicator to the etched surface and allows to dry for 30 seconds.
5. Apply Renew Universal / Total Etch Bonding adhesive and follow the instructions for bonding and restorative with composites or cements.

Warnings : Read the instructions thoroughly before use of this product. Hydrofluoric acid is a very powerful chemical and extreme care must be taken while handling this product.

Precautions :

1. Avoid direct exposure to sunlight.
2. Do not consume on take internally.

Contraindications: Do not use for etching of tooth surface (enamel or dentin).

Storage Condition : Store at temperature between 10°C(50°F) to 24°C(75°F). Close the cap carefully after use. Keep away from moisture. Shelf life not affected after opening if stored in a cool dry place.

Shelf Life : 2 years from the date of manufacturing.

Presentation :

Ref. 10012

Intro Pack : 1 x 1.2ml Syringes + 5 Applicator Tips.

Ref. 10013

Economy Pack : 1 x 10ml Plastic Bottle.

Disposal conditions : To be disposed off based on the prevailing rules of the country in which it is used.

Prevest DenPro Limited
Unit II, Export Promotion Industrial Park,
(EPIP), Bari Brahmna, Jammu-181133,
India. www.prestdenpro.com
Customer Care No. : + 91 8869074151
Email : customercare@prevestdenpro.com
Made in India

Obelis s.a
Boulevard General Wahis 53
1030 Brussels, Belgium
Tel: +32(0) 2 732 29 54 Fax: +32(0) 2 732 29 02
E-mail : mail@obelis.net

UDI Symbol for "UNIQUE DEVICE IDENTIFIER"

CC Symbol for "CORROSIVE"

MD Symbol for "MEDICAL DEVICE"

LOT Symbol for "BATCH CODE" **U** Symbol for "USE-BY DATE" **NON STERILE** Symbol for "NON STERILE" **KEEP DRY** Symbol for "KEEP DRY"

CAUTION Symbol for "CAUTION" **CONSULT INSTRUCTION FOR USE** Symbol for "CONSULT INSTRUCTION FOR USE" **REF** Symbol for "CATALOGUE NUMBER"

EC REP Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"

DO NOT USE IF PACKAGE IS DAMAGE Symbol for "DO NOT USE IF PACKAGE IS DAMAGE" **KEEP OUT OF THE REACH OF CHILDREN** Symbol for "KEEP OUT OF THE REACH OF CHILDREN" **MANUFACTURER** Symbol for "MANUFACTURER"

TEMPERATURE LIMIT Symbol for "TEMPERATURE LIMIT" **KEEP AWAY FROM SUNLIGHT** Symbol for "KEEP AWAY FROM SUNLIGHT" **DANGER** Symbol for "DANGER"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.