



PREVESTDenPro®

Reliable Disinfection

HYGIENE & DISINFECTION PRODUCTS

Medical Washers & Disinfectors



☒ FUNGICIDAL ☒ BACTERICIDAL ☒ TUBERCULOCIDAL ☒ VIRUCIDAL ☒ MYCOBACTERICIDAL ☒ YEASTICIDAL

www.prevesthygieneworld.com



The dynamic, forward looking, profession team of Prevest DenPro Limited Research & Development department has successfully developed a comprehensive range of revolutionary latest generation of broad spectrum, rapid action cleaning & disinfection solutions for dental and medical equipments. The products are recommended for use for cleaning and disinfection of all types of non invasive medical devices used in dental and medical operatories. The products have been tested for the bactericidal, virucidal, fungicidal, Sporicidal and yeasticidal efficacy. The products have been formulated with latest generation of patented blend of quaternary ammonium compounds, amines, surfactants and excipients which are highly effective against gram negative and gram positive bacteria, antibiotic resistant bacteria, human virus, veterinary viruses, fungi, tuberculosis. The products are free of phenol, aldehydes, halogens and chlorine. The product formulations are highly effective, easy to use and safe.

ViroGuard is ready to use fast acting and broad spectrum spray for quick cleaning and disinfecting of medical devices and surfaces after each session in medical and dental offices & clinics. It is specially designed to avoid cross contamination in between patients.

Key features of product :

- Bactericidal, fungicidal, tuberculocidal, mycobactericidal, virucidal
- Fast acting-Active in 30 seconds
- Does not affect the medical equipment
- Leaves no residue after drying
- Does not contain phenols, aldehydes, chlorine or EDTA

Indications for use:

It is used for cleaning, decontaminating and disinfecting surfaces such as trolleys, beds, chairs, trays, work benches, examination tables, dental units, handles, spittoons, other non invasive equipments and instruments.

Active Ingredients:

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2- Butoxyethanol Surfactant ,IPA , pH modifier and water.

Packaging:

REF. 16012, 500 ml spray bottle

REF. 16013, 1 Litre bottle

REF. 16014, 5 Litre bottle

Recommended usage:

Spray directly on the surfaces and objects to be disinfected and leave the surfaces wet for 30 seconds for effective disinfection and wipe with the clean lint free disposable cloth. (Note : Do not rinse)



ViroGaurd Foam

Alcohol Free Disinfectant Spray Foam For Medical Devices

ViroGuard Foam is an alcohol free formulation for rapid disinfection of alcohol sensitive non invasive medical devices and medical equipments having surfaces made of glass, rubber, acrylic plastic, and artificial leather upholstery.

Key features of product :

- Fast acting, ready-to-use alcohol free spray foam disinfectant
- Bactericidal, fungicidal, tuberculocidal, myco-bactericidal , virucidal
- Leaves no residue after drying
- Does not affect medical equipments
- Suitable for alcohol sensitive materials such as leatherite, plastic
- Does not contain phenol, aldehydes and chlorine
- Excellent wetting properties

Indications for use:

- It is an alcohol free foam spray used for disinfection of alcohol sensitive non invasive medical devices such as dental units, beds, stretchers, work benches etc.
- It can also be used for disinfection of non invasive medical equipments such as monitors, displays, keyboards and control panels etc.

Active Ingredients:

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2- Butoxyethanol Surfactant , D-Glucopyranose, oligomers, decyl octyl glycosides, Sodium lauroyl sarcosinate, pH modifier and water

Packaging:

REF. 16009, 500 ml spray bottle

REF. 16010, 1 Litre bottle

REF. 16011, 5 Litre bottle

Recommended usage:

Spray directly on the surfaces and objects to be disinfected and leave the surfaces wet for 30 seconds for effective disinfection and wipe with the clean lint free disposable cloth. For difficult surfaces such as keyboards, displays, monitors and control panels, use wet clean disposable cloth with spray foam and thoroughly wipe the surfaces.



ViroPerox is a powder concentrate for cleaning, cold sterilization and disinfection of all kinds of non invasive medical devices.

Key features of product :

- High efficiency concentrated formula for manual and ultrasonic cleaning.
- Bactericidal, fungicidal, tuberculocidal, virucidal and sporicidal.
- Exceptionally safe to use- Environment friendly
- Non- carcinogenic- Aldehyde and phenol free.

Indications for use:

- It is indicated for quick manual cleaning and disinfection of non-invasive medical devices.
- It can also be used in ultrasonic cleaners for cleaning and disinfection of instruments & equipments with increased risk of infection.

Active Ingredients:

Sodium Percarbonate, Tetraacetythylenediamine, Citric Acid Anhydrous, Disodium Metasilicate, Alcohol Ethoxylated, Sodium Carbonate, Trisodium Orthophosphate, Enzyme and pH modifier.

Packaging:

REF. 16007, 500g Jar

REF. 16008, 1 kg Jar

Recommended usage:

20g powder to be mixed with 1 litre of luke warm water. Soak the instruments for 10 minutes, remove, rinse with water and dry.



ViroZym is high efficiency concentrated liquid formula for manual and ultrasonic cleaning and disinfection of medical devices such as dental and medical instruments.

Key features of product :

- Effective with 1% dilution.
- Bactericidal, fungicidal and virucidal.
- Cleans, disinfects and deodorizes in one step.
- Removes stubborn deposits.
- Does not foam
- Does not contain phenols, aldehydes or chlorine

Indications for use:

- It is used for manual cleaning and disinfection of dental and medical devices such as rotary and surgical instruments respectively.
- It can also be used in ultrasonic cleaners for high level enzymatic cleaning and disinfection of dental and medical instruments.

Active Ingredients:

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine , Poly(oxy-1,2-ethanediyl),alpha-[2-(didecylmethyammonio)ethyl]- omega-hydroxy-,propanoate , Didecylidimethylammonium chloride, Isotridecanol, ethoxylated, 2-(2-Butoxyethanol), Subtilisin, thickening agents, emulsifiers and flavoring agents.

Packaging:

REF. 16005, 1 Litre bottle

REF. 16006, 5 Litre bottle

Recommended usage:

The concentrate is to be diluted with portable water to 1% concentration. Immerse instruments immediately after use in the solution and soak for 15 minutes. Rinse the instruments with water and dry after disinfection. Discard the used solution daily and prepare fresh solution daily. Do not mix with other cleaning agents or disinfectants.



ViroShield PDR

Powder Concentrate for Suction Unit Disinfection

ViroShield PDR is a powder concentrate for weekly hygiene and disinfection of dental suction systems of all types including their drainage pipes and spittoons.

Key features of product :

- Cost effective powder concentrate
- Bactericidal, fungicidal and virucidal
- Eradicates bad smell
- Does not foam
- Does not contain phenols, aldehydes or chlorine

Indications for use:

It is indicated for weekly hygiene and disinfection of dental suction systems as a part of their general disinfecting routine.

Active Ingredients:

Sodium Percarbonate, Tetraacetylenediamine, Citric Acid Anhydrous, Disodium Metasilicate, Alcohol Ethoxylated, Sodium Carbonate, Trisodium Orthophosphate, Enzyme.

Packaging:

REF. 16003, 500g Jar

REF. 16004, 1 kg Jar

Recommended usage:

20g powder to be mixed with 1 litre of luke warm water. Aspirate the solution through the high and low volume suction tubes. Leave for 30 minutes. Rinse thoroughly with clean water.



ViroShield is a specially formulated, foam free, liquid concentrate for daily cleaning, disinfection and deodorizing medical devices including suction devices and dental unit spittoons.

Key features of product :

- Cost effective-2% dilution
- Bactericidal, fungicidal and virucidal
- Eradicates bad smell
- Does not foam
- Does not contain phenols, aldehydes or chlorine

Indications for use:

It is indicated for daily hygiene and disinfection of dental suction systems as a part of their general disinfecting routine.

Active Ingredients:

Dimethyldioctyl ammonium chloride, N-(3-aminopropyl)-N-dodecylpropano-1,3-Diamine, Quaternary Ammonium Salts, Tetra Potassium Pyrophosphate, Isotridecanol Ethoxylated, 2-Aminoethanol, EDTA, Potassium Hydroxide, Dipropylene Glycol Methyl Ether and IPA.

Packaging:

REF. 16001, 1 Litre bottle

REF. 16002, 5 Litre bottle

Recommended usage:

The concentrate is to be diluted with water to 2% concentration. Pass the solution through the suction system once or twice a day. Leave for 15 minutes to take the effect. After disinfection, rinse the suction system with the portable water.



VIROPROTECT WIPES

Wipes for Rapid Disinfection
of Medical Devices

ViroProtect Wipes is a one-step, ready to use, disposable wipes for the rapid disinfection of the all kinds of wipeable surfaces of non-invasive medical devices and other non invasive instruments and equipments with increased risk of infection.

Key features of product :

- One-step-simply wipe to clean and disinfect surfaces of non-invasive medical devices.
- Revolutionary formula-impregnated with patented blend of antimicrobial compounds in alcoholic solution.
- Leaves no residue on the surfaces
- Fast acting- contact time 30 seconds
- Pleasant orange fragrance

Indications for use:

It is used for a rapid disinfection and cleaning of medical devices and surfaces like hand pieces, contra-angles, turbines, scaler inserts etc.

Active Ingredients:

N-(3-aminopropyl)-N-dodecylpropano-1,3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2- Butoxyethanol Surfactant , IPA, pH modifier, water

Packaging:

REF. 16015, 10 x wipes mono pack

Recommended usage:

Wipe the equipment with ViroProtect Wipe and massage thoroughly over the surface and provide scrubbing action to ensure good surface contact with the disinfectant. Leave the instruments/equipments to dry. Do not rinse with water and do not mix with any other disinfectant.



PREVESTDenPro[®]

Reliable Disinfection



Prevest DenPro Limited

Export Promotion Industrial Park,
Bari Brahmana, Jammu - 181133 (India)

+91 88990 74151 (Customer Care Number)
+91 88993 36171 (For International Market)
+91 98585 11041 (For Indian Market)

info@prevestdenpro.com
customercare@prevestdenpro.com

www.prevestdenpro.com