

# ViroPerox

**PREVESTDenPro®**  
THE FUTURE OF DENTISTRY

## Powder Concentrate for Medical Instruments Disinfection

**Product Description :** ViroPerox is a powder concentrate for cleaning, cold sterilization and disinfection of all kinds of non-invasive medical devices.

### Composition :

Sodium Percarbonate, Tetraacetylenediamine, Citric Acid Anhydrous, Disodium Metasilicate, Alcohol Ethoxylated, Sodium Carbonate, Trisodium Orthophosphate, Enzyme and pH modifier.

**Intended Purpose/Intended Use:** It is intended to be used for cleaning, decontaminating and disinfecting medical device and surfaces.

### Indications for use:

- ◆ It is indicated for quick manual cleaning and disinfection of non-invasive medical devices.
- ◆ It can also be used in ultrasonic cleaners for cleaning and disinfection of instruments & equipments with increased risk of infection.

**Instructions for Use :** 20g powder to be mixed with 1 liter of luke warm water. Soak the instruments for 10 minutes, remove, rinse with water and dry.

**Warnings :** Harmful if swallowed. Causes severe skin burns and eye damage. May cause damage to organs through prolonged or repeated exposure.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and continue rinsing. Get medical advice/attention if you feel unwell.

### Precautions :

- ◆ Advised for use on surfaces of medical devices by professionals only.
- ◆ Wear personal protective equipment such as gloves/protective clothing/eye protection/ masks for face protection. Wash your hands after handling.
- ◆ Do not mix with other products.
- ◆ If used to disinfect a large number of instruments heavily contaminated with blood, saliva or organic tissue, the solution should be changed more frequently.

### Contraindications : Avoid heat

**Storage Conditions :** Store at temperature between 10°C(50°F) to 24°C(75°F). Keep away from moisture and heat. Close the cap carefully after use. The shelf life of the product is valid only if left unopened and correctly stored in a cool dry place.

**Shelf life :** 3 years from the date of manufacturing.

### Presentation :

Ref: 16007

500g Jar

Ref: 16008

1kg Jar

**Disposal Conditions :** To be disposed off based on the prevailing rules of the country in which it is used.

### Prevest DenPro Limited

Unit II, Export Promotion Industrial Park,  
(EPIP), Bari Brahmna, Jammu-181133,  
India. [www.prevestdenpro.com](http://www.prevestdenpro.com)

Customer Care No. : + 91 8899074151

Email : [customercare@prevestdenpro.com](mailto:customercare@prevestdenpro.com)  
Made in India



Revision :	00
Date :	06-05-2023

UDI Symbol for "UNIQUE DEVICE IDENTIFIER"

LOT Symbol for "BATCH CODE" Symbol for "USE-BY DATE" Symbol for "NON-STERILE" Symbol for "KEEP DRY"

Symbol for "CAUTION" Symbol for "CONSULT INSTRUCTION FOR USE" Symbol for "CATALOGUE NUMBER"

Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"

Symbol for "DO NOT USE IF PACKAGE IS DAMAGED" Symbol for "KEEP OUT OF THE REACH OF CHILDREN" Symbol for "MANUFACTURER"

Symbol for "TEMPERATURE LIMIT" Symbol for "KEEP AWAY FROM SUNLIGHT" Symbol for "MEDICAL DEVICE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.