

ViroZym

PREVESTDenPro®
THE FUTURE OF DENTISTRY

Enzymatic Concentrate For Cleaning & Disinfection

Product Description : Virozym is a high efficiency concentrated liquid formula for manual and ultrasonic cleaning and disinfection of medical devices such as dental and medical instruments.

Composition :

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine , Poly(oxy-1,2-ethanediyl),alpha-[2-(didecylmethylammonio)ethyl]- omega-hydroxy-.propanoate , Didecylidimethylammonium chloride, Isotridecanol, ethoxylated, 2-(2- Butoxyethanol), Subtilisin, thickening agents, emulsifiers and flavoring agents.

Intended Purpose/Intended Use: It is intended to be used as a liquid formula for cleaning, decontaminating and disinfecting medical device and dental instruments.

Indications for use:

- ♦ It is used for manual cleaning and disinfection of dental and medical devices such as rotary and surgical instruments respectively.
- ♦ It can also be used in ultrasonic cleaners for high level enzymatic cleaning and disinfection of dental and medical instruments.

Instructions for Use : The concentrate is to be diluted with potable water to 1% concentration. Immerse instruments immediately after use in the solution and soak for 15 minutes. Rinse the instruments with water and dry after disinfection. Discard the used solution daily and prepare fresh solution daily. Do not mix with other cleaning agents or disinfectants

Warnings : Harmful if swallowed. Causes severe skin burns and eye damage. May cause damage to organs through prolonged or repeated exposure.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF ON SKIN(or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and continue rinsing. Get medical advice/attention if you feel unwell

Precautions :

- ♦ Wear personal protective equipment such as gloves, masks, goggles for protection. Wash your hands after handling.
- ♦ It is recommended to replace the solution, especially in the presence of high levels of organic contamination (blood , saliva)

Contraindications : Avoid heat

Storage Conditions : Store at temperature between 10°C (50°F) to 24°C (75°F). Keep away from moisture and heat. Close the cap carefully after use. The shelf life of the product is valid only if left unopened and correctly stored in a cool dry place.

Shelf life : 3 years from the date of manufacturing.

Presentation :

Ref: 16005

1 Litre Dosing Bottle

Ref: 16006

5 Litre Can

Disposal Conditions : To be disposed off based on the prevailing rules of the country in which it is used.

Prevest DenPro Limited

Unit II, Export Promotion Industrial Park,
(EPIP), Bari Brahmna, Jammu-181133,
India. www.prevestdenpro.com

Customer Care No. : + 91 8899074151

Email : customercare@prevestdenpro.com
Made in India



Revision :	00
Date :	25-04-2023

UDI Symbol for "UNIQUE DEVICE IDENTIFIER"

LOT Symbol for "BATCH CODE" Symbol for "USE-BY DATE" Symbol for "NON STERILE" Symbol for "KEEP DRY"

Symbol for "CAUTION" Symbol for "CONSULT INSTRUCTION FOR USE" Symbol for "CATALOGUE NUMBER"

Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"

Symbol for "DO NOT USE IF PACKAGE IS DAMAGED" Symbol for "KEEP OUT OF THE REACH OF CHILDREN" Symbol for "MANUFACTURER"

Symbol for "TEMPERATURE LIMIT" Symbol for "KEEP AWAY FROM SUNLIGHT" Symbol for "MEDICAL DEVICE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.