



PREVEST DENPRO LIMITED

Reg. Office : 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmana, Jammu-181133, India
Phones : +91 1923 222774, 222779, Fax : 222779
Web : www.prevestdenpro.com, E-mail : info@prevestdenpro.com

MATERIAL SAFETY DATA SHEET

SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

Product Name: Gingiva Shield VLC - Gingival Barrier
INTENDED USE: Dental
Company Name: Prevest Denpro Limited
Address: 38, Industrial Estate, Digiana, Jammu-180010 (India)
Tel: + 91-191-2430442
Fax No. +91-191-2430442
Email: info@prevestdenpro.com
Web: www.prevestdenpro.com
Emergency Contact No:- +91-191-2432548

SECTION 2: COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS#	EINECS	Classification	Content %
Diurethane Dimethacrylate	41137-60-4	255-239-5	Xi; R36, R37, R38	80.75

SECTION 3: HAZARDS IDENTIFICATION

HAZARD DESIGNATION:

"Xi"- Irritant

POTENTIAL HEALTH EFFECTS

Eyes: Causes eye irritation

Skin: May cause skin irritation.

SECTION 4: FIRST AID MEASURES

Skin Contact: Wash with soap and water.

Eye Contact: Hold eyelids apart and flush eyes with plenty of water for at least 15 minutes.

Ingestion: If swallowed, do **NOT** induce vomiting. Seek medical advice.

SECTION 5: FIRE-FIGHTING MEASURES

Fire Fighting Procedures:

General: Evacuate all personnel; use protective equipment for fire fighting. Use a self-contained breathing apparatus when the product is involved in fire.



PREVEST DENPRO LIMITED

Reg. Office : 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmana, Jammu-181133, India
Phones : +91 1923 222774, 222779, Fax : 222779
Web : www.prevestdenpro.com, E-mail : info@prevestdenpro.com

MATERIAL SAFETY DATA SHEET

Fire Fighting Equipment: Water spray, fog or mist, foam, dry chemical, carbon dioxide.

SECTION 6: ACCIDENTAL RELEASE MEASURES

Small Spill: Wipe up with cloth, soap and water.
Large Spill: Collect in tightly closed containers.

SECTION 7: HANDLING AND STORAGE

Handling: Keep away from sources of ignition and direct sunlight.
Storage: Store in a cool place in original container and protect from sunlight.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

PERSONAL PROTECTIVE EQUIPMENT:

Eyes and Face: Wear eye protection
Skin: S36/37: Wear suitable protective clothing.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Physical Form: Gel
Color: Blue green
Odor: Acrylic

SECTION 10: STABILITY AND REACTIVITY

Stability: Stable when stored and handled under recommended conditions
Polymerization: Polymerization occurs when exposed to visible light, ultraviolet light or extreme heat.

SECTION 11: TOXICOLOGICAL INFORMATION

Eye Effects: R36/38: Irritating to eyes and skin.

SECTION 12: ECOLOGICAL INFORMATION

Environmental Data: Not readily degradable.
Ecotoxicological Information: Water pollutant.



PREVEST DENPRO LIMITED

Reg. Office : 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmana, Jammu-181133, India
Phones : +91 1923 222774, 222779, Fax : 222779
Web : www.prevestdenpro.com, E-mail : info@prevestdenpro.com

MATERIAL SAFETY DATA SHEET

SECTION 13: DISPOSAL CONSIDERATIONS

Methods of Disposal: Dispose of in compliance with governmental regulation.
(EC 1975L0442-20/11/2003)

SECTION 14: TRANSPORT INFORMATION

Road and Rail (ADR/RID)

Label: Not classified

Air (ICAO/IATA)

Label: Not Classified

Vessel (IMO/IMDG)

Label: Not Classified

SECTION 15: REGULATORY INFORMATION

EUROPEAN COMMUNITY

ECC LABEL SYMBOL AND CLASSIFICATION:



"Xi" – Irritant

R36/37/38: Irritating to eyes, respiratory system and skin.

SECTION 16: OTHER INFORMATION

References: Not Available

Other Special Considerations: Not Available

Created: 27.04.2013

The information contained herein is based on data considered accurate. However, No warranty is expressed or implied regarding the accuracy of these data or the results to be obtained from the use thereof. In no event will the manufacture of the distributor be responsible for damages of any nature whatsoever resulting from the use of or reliance upon this information. General properties are to be regarded as guidelines and are not guaranteed for all samples