

Orafil-G™

Temporary Filling Material, Paste

Orafil-G™ is a self curing, ready to use, radiopaque cement for provisional filling of dental cavities. Orafil G is a zinc sulphate - zinc oxide based cement and is designed for short term temporary applications.

This product is also available with Sodium Fluoride as Orafil Plus



Indications

- Temporary restoration of cavities
- Temporary restoration between multiple endodontic visits
- Temporary restoration for inlay prepared cavity


Benefits

- Non eugenol formulation offers non irritating properties
- Resin addition gives high adhesive properties
- Quick and void free curing in moist environment
- Pulp protecting and bactericidal
- Slight expansion of filling material ensures well sealed margin
- Excellent adaptability
- Simple to apply with filling instrument


Chemistry

Zinc oxide, Zinc sulphate in EVA Resin/Glycol base with spearmint aroma


Application Technique




Endodontic Access



Placement of Material



Material Adaptation



Completed Temporary Restoration

Presentation

Orafil G™

Product Code : 11001
1 x Jar 40g, 1 x Instruction Manual
Product Code : 11002
6 xSticks 10g, 1 x Instruction Manual

Orafil Plus™

Product Code : 11003
1 x Jar 40g, 1 x Instruction Manual
Product Code : 11004
6 xSticks 10g, 1 x Instruction Manual

PREVEST DENPRO LIMITED

Unit I, 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmna, Jammu-181133, India

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com

Orafil-G™

Temporary Filling Material, Paste



Orafil-G™ is a self curing, ready to use, radiopaque cement for provisional filling of dental cavities. Orafil G is a zinc sulphate - zinc oxide based cement and is designed for short term temporary applications.

PREVEST DENPRO LIMITED

Unit I, 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmana, Jammu-181133, India

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com