

OraTemp NE™

Non-eugenol Temporary Luting Cement

OraTemp NE™ is a non eugenol temporary cement which permits placement and retention under function of permanent crowns, bridges, splints and their temporary counterparts. OraTemp NE is insoluble in oral fluids and retains a tight marginal seal to prevent leakage of these fluids under the restorations.



Indications

- Cementation of crowns, bridges and splints
- Retentions of temporary crowns and bridges
- Trial cementation of permanent restorations prior to final insertion

Benefits

- Excellent physical properties
- Easy to mix - easy to clean
- Low solubility in oral fluids
- Increased patient comfort
- Virtually odourless/tasteless

Chemistry

Ethoxy benzoic acid based, non eugenol temporary cement.
Contains: Ethoxy benzoic acid, Zinc oxide, Catalyst and Stabilizer

Application Technique



1 Dispense OraTemp NE



2 Spatulate to get a homogenous mix



3 Mix loaded on to a temporary restoration



4 Temporary crown cemented

Presentation

Product Code : 11003
1 x 15g Cartridge, Base & Catalyst,
1 x Spatula, 1 x Mixing Pad
1 x Instruction Manual

Product Code : 11004
1 x 15g Automix Syringe, Base & Catalyst
20 x Automix Tips
1 x Instruction Manual

PREVEST DENPRO LIMITED

Unit I, 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmana, Jammu-181133, India

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com

OraTemp NE™

Non-eugenol Temporary Luting Cement



OraTemp NE™ is a non eugenol temporary cement which permits placement and retention under function of permanent crowns, bridges, splints and their temporary counterparts. OraTemp NE is insoluble in oral fluids and retains a tight marginal seal to prevent leakage of these fluids under the restorations.

PREVEST DENPRO LIMITED

Unit I, 38, Industrial Estate, Digiana, Jammu-180010, India
Unit II, EPIP, Bari Brahmana, Jammu-181133, India

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com