

Pulp-X

Pulp Devitalizer

Pulp-X is a preparation for devitalization of the pulp. Pulp-X provides powerful pulp devitalizing properties due to paraformaldehyde.



Indications

- For Pulp devitalization.
- For residual devitalization after removal of non-vital pulp tissue.

Benefits

- Fast acting-pain releasing formula.


Chemistry

Formaldehyde Liquid, Parachlorophenol, Polyethylene Glycol, Amorphous Silica and Vinyl Siloxane.

Techniques



Removal of non vital pulp



Placement of Pulp-X



Seal hermetically with a temporary restoration

Presentation

- Ref. : 40021 (1 x 6g, Jar)
Ref. : 40022 (1 x 3g Syringe, 10 x Applicator Tips)
Ref. : 40022-1 (2 x 3g Syringe, 10 x Applicator Tips)

PREVEST DENPRO LIMITED

Export Promotion Industrial Park (EPIP),
Bari Brahmana, Jammu-181133, India.

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779
Toll Free No. : 1800-180-7172

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com
E-mail : customer@prevestdenpro.com

Pulp-X

Pulp Devitalizer



Pulp-X is a preparation for devitalization of the pulp. Pulp-X provides powerful pulp devitalizing properties due to paraformaldehyde.

PREVEST DENPRO LIMITED

Export Promotion Industrial Park (EPIP),
Bari Brahmana, Jammu-181133, India.

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779
Toll Free No. : 1800-180-7172

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com
E-mail : customer@prevestdenpro.com