

FUSION CRYSTA ORTHODONTIC ADHESIVE KIT™
Light Cure Orthodontic Kit

Crysta Orthodontic Adhesive Kit is a light curing composite kit used for bonding Metal & Ceramic brackets to tooth surface and for cementation of molar bands. The kit contains a revolutionary primer based on MDP technology. The novel technology ensures high enamel bond strengths. The combination of Crysta adhesive and primer provides reliable long lasting bonding in orthodontic cases.



Indications

- Bonding Metal & Ceramic brackets to tooth surface
- Cementation of molar bands.

Benefits

- Highly translucent adhesive designed to compliment traditional and LED curing lights
- Superior physical properties
- The primer contains the revolutionary 10-MDP
- MDP promotes physical and chemical adhesion
- High enamel bond strengths
- Excellent handling properties

Chemistry

FUSION CRYSTA ADHESIVE
Bisphenol a glycidyl methacrylate, Urethane dimethacrylate, Silane treated silica, Bisphenol a bis(2-hydroxyethyl ether) dimethacrylate, photopolymerizing agents.

FUSION CRYSTA ADHESIVE PRIMER
BisGMA, Triethylene Glycol Dimethacrylate, Camphorquinone, Dimethyl Amino Ethyl Methacrylate, Ethyl 4-dimethyl Amino Banzoate, 10-Methacryloyloxydecyl Dihydrogen Phosphate, Butylated Hydroxy Toluene.

Technique



Technical Data

- 1) Fusion Crysta Adhesive:**
Flexural strength: 120 Mpa
- 2) Fusion Crysta Primer:**
Bond strength : 29 Mpa

Presentation

Ref. : 10024

Kit Packing :

- ♦ 1 x 4g Syringe Fusion Crysta Adhesive
- ♦ 1 x 5ml Bottle Fusion Crysta Adhesive Primer
- ♦ 1 x 5ml Syringe Actino Gel
- ♦ 10 Applicator Tips
- ♦ 1 x Instruction Manual

PREVEST DENPRO LIMITED

Export Promotion Industrial Park (EPIP),
Bari Brahmna, Jammu-181133, India.

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779
Toll Free No. : 1800-180-7172

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com
E-mail : customercare@prevestdenpro.com

Fusion Crysta Orthodontic Adhesive Kit™ Light Cure Orthodontic Kit



Crysta Orthodontic Adhesive Kit is a light curing composite kit used for bonding metal & ceramic brackets to tooth surface and for cementation of molar bands. Novel MDP primer technology combined with the highly translucent adhesive, makes it a reliable combination for stress free orthodontic bonding

PREVEST DENPRO LIMITED

Export Promotion Industrial Park (EPIP),
Bari Brahmna, Jammu-181133, India.

Tel.: 91 191 2430442/ 2432548
Tel.: 91 1923 222774/222779
Toll Free No. : 1800-180-7172

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com
E-mail : customercare@prevestdenpro.com