


Florence®

In Office Whitening System

Florence is a chemically activated in office whitening gel containing 35% hydrogen peroxide. The new pH activated technology enables a fresh fast acting mix every time. This unique technology can be used to whiten teeth in a relatively short time. The system employs a novel Blockdam partially adhesive resin isolation technology to isolate the gingiva from the harmful effects of the peroxide ions. The whitening system does not require any light activation, making it the apt technology to be incorporated into every day dental practice.

Chemically activated in office
whitening system



Indications

- In office whitening for vital teeth.

Benefits

- No light activation.
- Relatively short application time.
- Novel isolation technology .
- Stabilized Hydrogen peroxide for better results and long shelf life.
- Easy application technique.

Chemistry

Hydrogen peroxide, alkali base.

Technical Data

1) PH : 9 - 10

Florence Bleaching Kit



Initial shade determination



Application of Liquid dam.
Light cure for 20 seconds



Application of Bleaching agent



Shade change after 20 minutes
of bleaching

Presentation

Ref. : 90002

Two Patient Kit :

- ◆ 4 x 0.3g Powder Pots
- ◆ 4 x 3ml Liquid Bottles
- ◆ 2 x 1.2ml Gingiva Shield VLC (Gingival Barrier Syringes)
- ◆ 1 x Mixing Pot
- ◆ 4 x Micro Brushes
- ◆ 4 x Dispensing Tips
- ◆ 1 x Instruction Manual

Ref. : 90001

Single Patient Kit :

- ◆ 2 x 0.3g Powder Pots
- ◆ 2 x 3ml Liquid Bottles
- ◆ 1 x 1.2 Gingiva Shield VLC (Gingival barrier Syringes)
- ◆ 2 x Micro Brushes
- ◆ 2 x Dispensing Tips
- ◆ 1 x Instruction Manual

PREVEST DENPRO LIMITED


Unit II, Export Promotion Industrial Park (EPIP),
Bari Brahmana, Jammu-181133, India.

Tel.: + 91 191 350 6858, 60, 61
Tel.: + 91 8899074151

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com

Florence®

In Office Whitening System



Chemically activated in office
whitening system

Florence is a chemically activated in office whitening gel containing 35% hydrogen peroxide. The new pH activated technology enables a fresh fast acting mix every time. This unique technology can be used to whiten teeth in a relatively short time. The system employs a novel Blockdam partially adhesive resin isolation technology to isolate the gingiva from the harmful effects of the peroxide ions. The whitening system does not require any light activation, making it the apt technology to be incorporated into every day dental practice.

PREVEST DENPRO LIMITED

Unit II, Export Promotion Industrial Park (EPIP),
Bari Brahmana, Jammu-181133, India.

Tel.: +91 191 350 6858, 60, 61
Tel.: +91 8899074151

E-mail : info@prevestdenpro.com
Website : www.prevestdenpro.com